

JONARIN

The Chemical People

Since 1976

An ISO 9001-2000 Certified Company

JONARIN'S SPANDAGE



WE BRING THE BEST FOR YOU FROM THE WORLD AROUND



JONARIN'S SPANDAGE

WHAT IS IT?

An open, webbed tubular shaped retainer stretch net composed of cotton, polyester fiber and high elastic, used for secure wound dressing. Tapes and adhesives are not required thus preventing allergy / irritation. The flexible net design of JONARIN's SPANDAGE allows free breathing and movement.



USE

- Efficient wound management easy examination and dressing change.
- Especially effective in burn care, surgeries, post-op dressing, donor wound graft sites
- Daily dressing changes.
- Retainer net holds hot and cold packs, iv tubes, medical devices, pads used for incontinence, diabetic ulcer & chronic wounds.
- SPECIFIC INDICATION - Patients with fragile skin and allergic to adhesive tapes.

MECHANISM OF ACTION

The gentle pressure provided by our SPANDAGE ensure proper microcirculation enabling FASTER wound healing than other bandages. The open web structure ensure optimum ventilation thereby ELIMINATING itching and skin discoloration.

WHY JONARIN'S SPANDAGE & NOT BANDAGE?



JONARIN'S SPANDAGE stretches Triple its original size ie, less material consumption!! (Cost effective than any conventional dressing material). All thread joints are knitted you can cut opening to any size without the risk of becoming damaged. Spandage does not obstruct circulation so does not cause skin irritation.

Thus, ENSURE 100% COMFORT WITH 0% EFFORT.


FEATURES & BENEFITS.

1. Washable, reusable, flexible and aerated
2. Fit easily over any body part (lower abdominal, groin, joints) & also over bulky dressing.
3. No use of adhesive tapes other fastening (no pain from tape removal, no tape excoriation problem & no fear of allergies. Fits well on hairy skin)
4. Allows easy wound inspection and quick application for wounds which need frequent dressing.
5. Comfortable and secure, absolute freedom of movement.
6. Promote micro blood circulation.
7. Non irritant & non allergic (ideal for all skin types skin).
8. Can be used on SWOLLEN WOUNDS due to its stretching ability.
9. Cuts easily & stretches to various body parts and shapes.
10. Used as dressing retention.

WARNING: Do not use 2 layers of dressing to any appendage.

Head - Full Cap

- Use size 7.
- Start with a 10cm or 12cm piece of Stretch Net.
- Seal one end of the Stretch Net by tying a knot 4cm from the open end, make a 1cm cut.




1. Place stretch Net over the head with the hole facing forward. Stretch the hole over the chin.



2. Complete dressing.

CHEST

- Use size 7 - 8 depending on patient size.
- Cut a 24cm piece of stretch Net. 4cm from one end, make a $\frac{1}{2}$ cm cut on both sides.



1. Starting with the end farthest from cuts, gather Stretch Net and slip over head. Make sure to keep cuts to the sides.




2. Continue to Stretch the Net over the shoulders and down the trunk.



3. Pass arms through holes. Complete dressing.

SHOULDER

- Use size 7 - 8
- Cut a 24cm piece of Stretch Net. 2cm from one end, make a 6cm slit through one side of the tube.




1. Starting with the end closest to the slit, slide stretch Net onto arm, with cut facing up.



2. Pull the cut loop over the head let rest under the opposite axilla. Complete dressing.

AXILLA OR BREAST

- Use size 7 - 8 depending on the size of the patient.
- Cut a 24cm piece of Stretch Net. 14cm from one end, make a 3 cm slit long ways in one side of the tube.



1. Starting with the end farthest from the slit, pass the Stretch Net over the arm up to the shoulder. Make sure to keep the cut facing up. Pass the cut loop over the head.




2. Continue stretching the cut over the opposite shoulder and arm allowing it to rest under the opposite axilla.

Complete dressing.

GENITAL & ANAL

- Use size 8 depending on patient size.
- Cut a 24cm piece of Stretch Net. In the middle of this piece make a 1cm V cut.



1. Beginning where the cut was made, insert both hands into one half of the tubing and stretch it sideways.



2. Insert both feet into the stretched half and pull it up over the waist. The second half should be in the back hanging like a tail.



3. Pull the tail between the legs and pull up over the genital area in a spreading motion.



4. Secure ends by tying them off in the front for a completed dressing.



5. Back view.

GROIN & BUTTOCKS

- Use size 7 - 8
- Cut a 14cm piece of Stretch Net. About 4cm from one end, make a 1/2cm cut.



1. Put your hands into the tube in the end farthest from the cut. Allow your left hand to go through the cut and your right hand to go through the longer tube.



2. Now stretch to form a panty.



3. Slip patient into the panty.

HIP OR THIGH

- Use size 5 or 7
- Cut a 24cm piece of Stretch Net. 4cm from one end, make a 1cm cut into the tube.



1. Starting with the end closest to cut, insert left hand into long portion of the tube and right hand into the tube down to the 1cm cut. Now stretch apart to create a waist band.



2. Slip affected leg into Stretch Net.




3. Stretch the waist band open and stretch it over the unaffected leg.



4. Slip Stretch Net up to waist for complete dressing

FOOT

- Use size 3 or 5
- Start with a 12cm piece of relaxed Stretch Net.



1. Slip Stretch Net onto the toes. Twist 2 times and then slip remaining Stretch Net over entire foot.




2. Clip threads to release pressure on heel.



3. Complete dressing.

ELBOW, FOREARM, OR KNEE

- Use size 3 or 5
- Start with a 12cm piece of relaxed Stretch Net.




1. Insert both hands into tube and stretch over elbow to complete dressing.



2. Insert both hands into tube and stretch over knee to complete dressing.

HAND

- Use size 3
- Start with a 12cm piece of relaxed Stretch Net.



1. Slip Stretch Net over hand with the palm facing upward. Leave 3cm extension beyond fingers.




2. At the tip of each finger cut 2 or 3 threads to release fingers.



3. Pull the extension down and tuck it under. Complete dressing.

FINGER

- Use size 5
- Cut a 12cm piece of relaxed Stretch Net. 2cm from one end, make a 1 1/2cm slit on one side.



1. Starting with end farthest from slit, slip Stretch Net on finger with slit facing down.



2. Twist the remaining Stretch Net one complete turn until the slit faces down again. Gather the remaining Stretch Net up to the finger tip.



3. Pull Stretch Net back over finger, slipping hand through slit, and anchor to wrist. complete dressing.



JONARIN PIGMENTS PVT. LTD

1st Floor, Vettikkat Tower, National Highway
Palarivattom, Kochi - 682 025

Tel: Off: 0484 2345556, 2346662, 2346044

Fax: 0484 2331192, Grams: jonarin

Email: jonarin@reliancemail.net, jonarin@vsnl.com

Website: www.jonarinjaliz.com